




ULTRASONIC CLEANER BULC-919

ULTRASONIC CLEANER BULC-919

ULTRASONIC CLEANER

Ultrasonic cleaner is a non-polluting, non-damaging, multi-purpose cleaning equipment. It not only has a powerful cleaning effect on objects with complex shapes, but also accelerates dissolution and emulsification. Ultrasonic cleaner has the unparalleled advantages of other physical cleaning or chemical cleaning, and is widely used in service industry, electronics industry, laboratory, machinery industry, cemented carbide industry, chemical industry and other fields.



1. Standard accessories: 1pc of soundproof cover, 1pc of basket, 1pc of PVC pipe.
2. The shell and inner cavity are made of 304 stainless steel, which is anti-corrosion.
3. LED display

SPECIFICATIONS

Model	BULC-919
Capacity (L)	22.5
Ultrasonic Frequency	Single frequency, 40kHz
Time Range/Heating Temperature	1~99min/RT~80°C
Ultrasonic Power	480W
Heating Power	500W
Power Supply	220V, 50Hz
Tank Size(mm)	500*300*150
External Size(mm)	530*330*310
Packing Size(mm)	640*415*380
Gross Weight(kg)	14
Alt Name	Ultrasonic Cleaner



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com